



MISSOURI DEPARTMENT OF HEALTH
STATE PUBLIC HEALTH LABORATORY
BAC VERIFIER MAINTENANCE REPORT

Complete this report in duplicate at the time of the regular monthly preventive maintenance check, and whenever instrument is repaired. Send copy to Department of Health; retain original in department file.

BAC VERIFIER SN	DATE OF INSPECTION
LOCATION OF INSTRUMENT (STREET AND CITY)	TIME OF INSPECTION

CHECKLIST: Place a check (✓) to the left of each item if found to be satisfactory or if operating within established limits. (Write in observed values where determined.) Unchecked items must be corrected before using instrument.

☐ DIAGNOSTIC CHECK (PRINTOUT ATTACHED)

☐ COMPUTER

☐ DETECTOR

☐ PROGRAM

☐ FILTERS

☐ HEATERS SAMPLE CHAMBER

☐ ACETONE SWITCH

☐ SET _____ °C

☐ QUARTZ STANDARD

☐ ACTUAL _____ °C

☐ CALIBRATION

☐ PUMP HIGH SPEED

☐ PRINTER

☐ INDICATOR LIGHTS

☐ TIME AND DATE

☐ SIMULATOR TEMPERATURE (34 °C ± 0.2°C)

☐ CALIBRATION CHECK -

Run three tests using a standard solution. All three tests must be within ± 5% of the standard value and must have a spread of .005 or less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)

☐ 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE

☐ 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE

(ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT)

TEST 1

TEST 2

TEST 3

☐ PERFORM R.F.I. TEST (PRINTOUT ATTACHED)

☐ NUMBER OF REFUSALS, SINCE LAST MAINTENANCE REPORT, AND NUMBER OF BREATH TESTS IN EACH RANGE AS FOLLOWS: **(DO NOT INCLUDE SIMULATOR TESTS)**

REFUSALS

(0-.04)

(.05-.09)

(.10-.14)

(.15-.19)

(Over .19)

List any new parts and describe any alteration or modification that was made to restore the instrument to operate satisfactorily and within established limits (use other side if necessary)

INSPECTING OFFICER

SIGNATURE



PRINT NAME

TYPE II PERMIT NUMBER/EXPIRATION DATE

TELEPHONE NUMBER